



Catatan Hasil Kalibrasi Internal
Internal Calibration Note
Moisture Analyzer Balance

No. : F-PM-01-49

Rev : 00

Berlaku : 03 September 2014

No. Fasilitas / Facility No. :
 No Kode Kalibrasi / Cal. Code :
 Merek / Brand :
 Tipe / Type :
 No. seri / Serial Number :
 Lokasi & Bidang / Location & Dept :
 No. Protap Cara Kalibrasi / SOP Calibration :
 Suhu Ruangan / Room Temp. :
 Kelembaban Ruangan / Room Humidity :

HASIL PENGUJIAN BOBOT / Weighing test

No.	Standar Bobot / Weight (g)	Penunjukan / Reading (g)	Kesalahan / Error (g)	Toleransi / Tolerance (g)	Keterangan / Note
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					

HASIL PENGUJIAN SUHU / Temperature test**Setting Temperature :**

Waktu / Time (min.)	1	2	3	4	5	6	7	8	9	10
Suhu Alat / Temp. Display (°C)										
Suhu Termokopel / Temp. Thermocouple (°C)										

Catatan / Note :

1. Anak Timbangan yang digunakan / Weight set Calibrator :
2. Kalibrator suhu yang digunakan / Temp. Calibrator :
3. Tgl Kalibrasi kalibrator / Cal. Date of Calibrator :
4. Tgl Rekalibrasi kalibrator / Recal. Date of Calibrator :
5. Metoda pengujian / SOP :
6. Tanggal pengujian / Cal date :
7. Disarankan kalibrasi ulang / Recal. date :

Bekasi,

.....
Ttd Petugas / Sign......
Sign. Of Asst. Manager