

	Catatan Hasil Kalibrasi Internal <i>Internal Calibration Note</i> Granulate Flow Tester	No. : F-PM-01-63
		Rev. : 00
		Berlaku/ Date : 3 September 2014

No. Fasilitas / <i>Facility No.</i>	:	Lokasi / <i>Location</i>	:
No. Kode Kalibrasi / <i>Calibration code</i>	:	Bidang / <i>Department</i>	:
Merek / <i>Brand</i>	:	Suhu ruangan / <i>Room temp.</i>	:
Tipe / <i>Type</i>	:	Kelembaban udara / <i>RH</i>	:
No seri / <i>Serial Number</i>	:	No. Protap Cara Kalibrasi / <i>SOP</i>	:
Kapasitas / <i>Capacity</i>	:	Nama Petugas / <i>Name</i>	:

Kalibrator yang digunakan <i>Calibrator</i>	Tanggal Kalibrasi Kalibrator <i>Cal. date of Calibrator</i>	Tanggal Rekalibrasi Kalibrator <i>Recal. date of Calibrator</i>

1. Kalibrasi pengukur waktu (satuan gram/detik) / Calibration of Timer (g/sec) Tanggal/ Date :

No	Bobot/ Weight (gram)	Pembacaan waktu/ Time Reading (detik/second)		Selish/ Difference	Syarat/ Requirement	Kesimpulan/ Conclusion
		Alat/ Instrument	Stopwatch/ Reference			
1					±1 detik/ second	
2						
3						
No	Bobot/ Weight (gram)	Pembacaan waktu/ Time Reading (detik/ second)		Selish/ Difference	Syarat/ Requirement	Kesimpulan/ Conclusion
		Alat/ Instrument	Stopwatch/ Reference			
1					±1 detik/ second	
2						
3						
Paraf Petugas/ <i>Signature</i>						
Paraf Asman/ <i>Sign. of Asst. Manager</i>						

	Catatan Hasil Kalibrasi Internal <i>Internal Calibration Note</i> Granulate Flow Tester	No. : F-PM-01-63
		Rev. : 00
		Berlaku/ Date : 3 September 2014

2. Kalibrasi pengukur waktu (satuan ml/detik) / Calibration of Timer (ml/sec)
Tanggal/ Date :

No	Volume (ml)	Pembacaan waktu/ Time Reading (detik/second)		Selish/ Difference	Syarat/ Requirement	Kesimpulan/ Conclusion
		Alat/ Instrument	Stopwatch/ Reference			
1					±1 detik/ second	
2						
3						
No	Volume (ml)	Pembacaan waktu/ Time Reading (detik/second)		Selish/ Difference	Syarat/ Requirement	Kesimpulan/ Conclusion
		Alat/ Instrument	Stopwatch/ Reference			
1					±1 detik/ second	
2						
3						
Paraf Petugas/ Signature						
Paraf Asman/ Sign. of Asst. Manager						